

RC-Base AbutmentTM

with MEGA¹GEN

AnyRidge made
one more
Big innovation!



Another innovation in AnyRidge prosthetics for anterior aesthetics

Stable crown fixation
structure

- Cement contact area ↑
- Wall thickness ↑

0.6mm Butt joint design

Prevents chipping of margins

Innovative
DP Screw™ (Double-Pitch Screw)

Unique design with two pitches

- Minimizes sinking
- Maintains occlusion relationship of prosthesis

Optimal angle option
for aesthetic zone

→ 0°/ 11°/ 22°

- No screw hole in aesthetic area

Ideal cementation surface

Improved cementation with ideal
surface geometry

Anti-rotation structure

Prevents any rotation movement
with crown

Biologic S-line

Ensures a better peri-implant
biotype & better emergence profile

Ideal angled screw channel

No screw hole in aesthetic area ensures best aesthetic result.

Perfect surgery and ideal implant position but...



RC-Base abutment™

Conventional abutment

· The angled tool channel allows the
screw hole to be removed from the
aesthetic area.



· The access channel for a screw-
retained prosthesis can interfere
with the aesthetics if the implant
is angled too far labially. In such
cases, a cemented type of retention
is the usual choice.



RC-Base abutment™ Conventional abutment


Innovative DP Screw™ (Double-Pitch Screw)

Structured to minimize sinking, which is a key disadvantage of all
internal fixtures




Stronger cementation force with zirconia custom

Maximized contact surface area



- The surface roughness of the RC-Base abutment™ is created by sandblasting.
- This maximizes the contact surface area with cement.
- Prevents cement overflow close to the screw channel and near the platform to the basal face.


Wider zirconia wall thickness due to better & easier screw driver access



- Better wall thickness is secured via angled tool channel.
- Abutment is provided with screw already attached, so no space needed for screw head.



No side effects from residual cement



① Final crown is cemented outside oral cavity ② All excess cement is removed from final prosthesis

RC-Base can solve your worries about ...


- If you want to make a zirconia custom abutment with CAM equipment
- If you are concerned about a screw hole in the aesthetic area
- If you are worried that the contour formation will not work properly
- If you are concerned about the strength of your zirconia custom abutment
- If you are concerned about screw loosening due to excessive mastication
- If there is a risk of inflammation due to residual cement in the case of deep placement

➡ Abutment options

RC-Base Abutment (Straight)

- DP Screw(1-MTO-D2T5) included
- Use with Screw Driver(MTO-DIT5-FW)
- Various cuff heights (2/ 3/ 4mm)
- Recommended tightening torque: 25Ncm

Diameter	Angle	Post Height (mm)	Cuff Height (Labial/Lingual) (mm)	Ref.C
Ø4.5	0°	4.7	2 (1.1)	D-MA35-C0G2-AS
			3 (1.9)	D-MA35-C0G3-AS
			4 (2.8)	D-MA35-C0G4-AS



RC-Base Abutment (11°)

- DP Screw(1-MTO-D2T5) included
- Use with Screw Driver(MTO-DIT5-FW).
- Various cuff heights (2/ 3/ 4mm)
- Recommended tightening torque: 25Ncm

Diameter	Angle	Post Height (mm)	Cuff Height (Labial/Lingual) (mm)	Ref.C
Ø4.5	11°	4.7	2 (1.9/1.1)	D-MA35-C1G2-AS
			3 (2.7/1.9)	D-MA35-C1G3-AS
			4 (3.6/2.7)	D-MA35-C1G4-AS



RC-Base Abutment (22°)

- DP Screw(1-MTO-D2T5) included
- Use with Screw Driver(MTO-DIT5-FW)
- Various cuff heights (2/ 3/ 4mm)
- Recommended tightening torque: 15Ncm

Diameter	Angle	Post Height (mm)	Cuff Height (Labial/Lingual) (mm)	Ref.C
Ø4.5	22°	4.7	2 (2.8/1.1)	D-MA35-C2G2-AS
			3 (3.5/1.9)	D-MA35-C2G3-AS
			4 (4.4/2.7)	D-MA35-C2G4-AS



DP Screw

- Double-Pitch Screw

Torx	Ref.C
T5	1-MTO-D2T5



Screw Driver

- Fracture torque : 35Ncm

Torx	Length (mm)	Ref.C
T5	30	MTO-DIT5-FW-BOX



Try-In Tool

• Used to move RC-Base abutment to fixture

Type	Length (mm)	Ref.C
M1.8	30	MTO-DT-BOX



Right Angle Adapter

• Connect to Screw Driver(MTO-DIT5-FW) for use with Torque Wrench

Type	Ref.C
Handpiece type	TTA100






➡ Components of Trial Intro Set

Ref.C
MA35-Set-16



➤ RC-Base Clinical Workflow




Fixture Positioning Guide
For the angled-type RC-Base abutment™, the post is inclined toward the edge of the connection, so the edge direction should be in the labial (Lingual) direction when installing the fixture.
※ The hex structure of the handpiece connector allows the hex direction of the fixture connector to be adjusted according to the implant placement.



Oral scan using scan abutment
Insert the scan abutment in accordance with the hex direction of the inner surface of the fixture, and connect the abutment screw. Check that the scan abutment is completely connected to the inner surface of the fixture, then perform oral scan.


» RC-Base Clinical Workflow

Optional flow





Try-In Tool

After separating the scan abutment, use the Try-In Tool to check the cuff height and angulation of the RC-Base abutment™. Use the Try-In Tool for easy pick-up and connection. Connect the screw tip of the Try-In Tool to the innerthread of the top part of the RC-Base abutment in a clockwise direction.




RC-Base Try-In

Check the cuff height and the position of the screw hole by connecting the RC-Base abutment™ attached to the Try-In Tool to the implanted fixture. The position of the screw hole can be adjusted in 6 directions based on the hex structure of the RC-Base abutment. Try to position the direction of the screw hole on the lingual side. For the accurate prosthetic work, take a photo of the screw hole direction at this time, or use the anti-rotation structure of the post to remember the direction.




CAD Order Sheet and Library Selection

After running the CAD (exocad, 3Shape) program, create an order sheet for the design work. While the order sheet format will differ according to the program (exocad, 3Shape), a screw-retained crown is the normal selection. Select the RC-Base library that matches the cuff height and angulation determined at the Try-In. Download and install the exocad or 3Shape library from the MedTEOR website or R2GATE website. CAD works well with the 2019 version or a later version.




Matching

After importing the scan data, start the normal matching process. As matching references, use the flat surface of the Scan Abutment and flat surface of the Scan Abutment in the scan data. After matching, check the degree of precision using the recommended function in respective CAD software.





Confirmation of Library, Direction, Cuff Height, and Diameter of RC-Base Abutment™

After matching the library, check whether the post direction of the replaced RC-Base matches the direction recorded during the Try-In. Also, check whether the cuff height and diameter of the selected library are appropriate. If not, return to the order sheet screen and select the correct library.




Prosthetic Design

After setting the margins according to the general rule, complete the design. Check that the screw hole is properly set in an aesthetic position and there are no weak areas in the crown design. In particular, make sure the thickness of the prosthesis is adequate near the lingual screw hole.



Extraction of STL

Extract or save the completed crown design in an STL format for importing into the CAM software. Set the desired processing area and perform milling and tool simulation according to the operation of the CAM software






Zirconia Milling




Import the extracted prosthesis STL file into the CAM software and set the desired area needed for processing. Perform milling and tool simulation according to the operation of the CAM software. Save the NC file in the CAM software and input the file into the milling equipment. Set the designated zirconia block and start machining.

>> RC-Base Clinical-work flow

RC-Base Abutment™

Conventional Ti-base Abutment



For the same anterior case, the screw hole designed for the RC-Base abutment™ is located on the lingual side, which results in an aesthetic and clean look on the labial surface.



For the same anterior case, the screw hole designed for a regular Ti-base abutment is located on the labial side, which makes it difficult to ensure a good aesthetic result.

After milling





Baking



Cementing





Picking up screw-retained prosthesis

Sintering Zirconia
 Sinter the finished zirconia crown according to the recommendations for the zirconia block. Post-process as normal and dye, if required.

Cement Bonding - Prosthesis & RC-Base
 Before cementing, combine the prosthesis and RC-Base and check the gap. Rework if the prosthesis is not properly coupled or it rotates due to excessive clearance.
 According to the cement manufacturer's instructions, apply the recommended amount of cement to the post of the RC-Base, then slowly attach it to the prosthesis. Remove any excess cement overflow using a sculpting tool. To complete, light-cure or self-cure according to the cement manufacturer's instructions.


Assembly using Screw Driver
 To fasten the prosthesis, connect the Right Angle Adapter to the Screw Driver or just use a Right Angle Driver.

Picking up the Prosthesis
 Pick up the prosthesis by connecting the Screw Driver to the DP Screw™.

10

11

Rotate & pull up screw until it reaches top of internal screw channel of RC-Base



Initialize Position of DP Screw™

To fasten the RC-Base to the fixture, rotate and pull up the DP Screw™ in a counter clockwise direction until the stopper of the DP Screw™ touches the lower part of the inner thread of the RC-Base.

※ Although the position of the DP Screw™ of the RC-Base is set at the top, it can move during the delivery stage or for other reasons. So, be sure to check the position of the DP Screw™ to achieve a complete connection with the fixture.

Transport of Prosthesis

Move the prosthesis using the Screw Driver.



Delivery of Prosthesis

According to the connection direction of the fixture and the adjacent teeth, turn the DP Screw™ in a clockwise direction to fasten the prosthesis to the fixture.

※ If it feels too tight initially, this is probably because the DP Screw™ is in contact with the internal thread of the RC-Base when fastening the RC-Base to the fixture. No problem, just rotate in a clockwise direction to release the stress.

※ Recommended torque: 25N(0°/11°)/ 15N(22°)



Delivery of prosthesis.
Recommended tightening
torque is 25Ncm.

→ Screw hole in lingual direction

Final prosthesis is complete, and since the screw hole is on the lingual side, this creates a better aesthetic result.



Rev. 01