

Custom abutment milled within 13 minutes


TiGEN ABUTMENT™ [cuff type]

- Minimizes patient pain when fitting temporary abutment, as size is matched with healing abutment
- Pre-milled cuff creates perfect fit with emergence profile, making it easy, even for beginners
- Optimized chairside solution with dramatically reduced milling time

Cuff-type TiGen Abutment™ for easy & fast processing of custom abutment

1. Pre-milled cuff sizes that match healing abutment to create perfect fit with emergence profile

- Various cuff sizes for every gingiva height



2. Pre-milled cuff reduces milling time + precision is increased with reverse jig milling

- >60% reduction in milling time when compared with conventional products!
- NO post milling, allowing reverse jig milling to occlusal surface within 8 minutes!



3. Exclusive library, making design easy, even for beginners

- Optimized for immediate loading
- Easy implant prosthesis design by cuff type

